

# In vivo 7 T diffusion MRI: High spatio-angular-temporal resolution pursuit

Zhengguo Tan <sup>1</sup>, Patrick A Liebig <sup>2</sup>, Robin M Heidemann <sup>2</sup>, Frederik B Laun <sup>3</sup>,  
Florian Knoll <sup>1</sup>

<sup>1</sup>Artificial Intelligence in Biomedical Engineering, University of Erlangen-Nuremberg, Erlangen, Germany

<sup>2</sup>Siemens AG, Erlangen, Germany


<sup>3</sup>Institute of Radiology, University Hospital Erlangen, Erlangen, Germany

07 June 2023

## Motivation: In vivo Diffusion MRI at 7 T

# MUSE<sup>1</sup> faces challenges in ultra-high field

- ✓ MUSE: self-gated multi-shot interleaved EPI (iEPI) DWI reconstruction




- \* requires fully-sampled in-plane acquisition

<sup>1</sup>Chen NK, et al. A robust multi-shot scan strategy for high-resolution diffusion weighted MRI enabled by multiplexed sensitivity-encoding (MUSE). *NeuroImage* (2013). doi: [10.1016/j.neuroimage.2013.01.038](https://doi.org/10.1016/j.neuroimage.2013.01.038)

# MUSE<sup>1</sup> faces challenges in ultra-high field

- ✓ MUSE: self-gated multi-shot interleaved EPI (iEPI) DWI reconstruction




- \* requires fully-sampled in-plane acquisition
- \* only limited number of shots have been reported

---

<sup>1</sup>Chen NK, et al. A robust multi-shot scan strategy for high-resolution diffusion weighted MRI enabled by multiplexed sensitivity-encoding (MUSE). *NeuroImage* (2013). doi: [10.1016/j.neuroimage.2013.01.038](https://doi.org/10.1016/j.neuroimage.2013.01.038)

# MUSE<sup>1</sup> faces challenges in ultra-high field

- ✓ MUSE: self-gated multi-shot interleaved EPI (iEPI) DWI reconstruction



- \* requires fully-sampled in-plane acquisition
- \* only limited number of shots have been reported
- \* increases specific absorption rate (SAR) at ultra-high field

<sup>1</sup>Chen NK, et al. A robust multi-shot scan strategy for high-resolution diffusion weighted MRI enabled by multiplexed sensitivity-encoding (MUSE). *NeuroImage* (2013). doi: [10.1016/j.neuroimage.2013.01.038](https://doi.org/10.1016/j.neuroimage.2013.01.038)


# Undersampled iEPI Acquisition and Joint $k$ - $q$ -Space Reconstruction

This work presents preliminary results on

1. undersampled iEPI acquisition to reduce scan time and to reduce SAR
2. joint  $k$ - $q$ -space reconstruction to boost SNR


## Methods

# Undersampled iEPI with Shift Encoding



**(A)** Given the acceleration along  $k_y$  as  $R \rightarrow (R \times N_{\text{shot}})$ -fold undersampling per shot


# Undersampled iEPI with Shift Encoding



- (A) Given the acceleration along  $k_y$  as  $R \rightarrow (R \times N_{\text{shot}})$ -fold undersampling per shot
- (B)  $k_y$  shift encoding along diffusion encoding  $\rightarrow$  complementary  $k$ - $q$ -space sampling

# In Vivo Acquisition

- ▶ three-shell diffusion sampling
- ▶ 1.4 mm isotropic resolution and 68 slices
- ▶  $k_y$  undersampling  $R = 3$ , partial Fourier 6/8
- ▶ 1- and 2-shot scans



## Joint $k$ - $q$ -Space Reconstruction

Shot phase variation informed joint reconstruction:

$$\operatorname{argmin}_{\mathbf{x}_q} \left\| \mathbf{y}_{c,q,s} - \mathbf{P}_{q,s} \mathbf{F} \mathbf{S}_c \{ e^{i\Phi_{q,s}} \mathbf{x}_q \} \right\|_2^2 + \lambda \mathbf{R}(\mathbf{x}_q) \quad (1)$$

$\mathbf{x}$  : diffusion weight images,

$\Phi$  : shot-to-shot phase variation,

$\mathbf{S}$  : coil sensitivity maps,

$\mathbf{F}$  : 2D Fourier transform,

$\mathbf{P}$  : sampling mask,

$\mathbf{y}$  : measured raw data.

Indices:  $q$  - diffusion encoding;  $s$  - shot;  $c$  - coil

# Overlapping Locally Low Rank (LLR) Regularization


1. Stride set as 1, leading to overlapping local patches
2. Use local patches to create spatial-contrast (diffusion) matrices
3. singular value thresholding (SVT)<sup>a</sup> on these matrices
4. Backpass (adjoint) the SVT matrices back to full image matrices

---



<sup>a</sup>Cai JF, Candes EJ, Shen Z. A singular value thresholding algorithm for matrix completion. *SIAM J Optim* (2010).  
doi: [10.1137/080738970](https://doi.org/10.1137/080738970)

## Results

# 1-Shot EPI: LLR Regularization Strength




## 1-Shot EPI: Overlapping LLR




## 2-Shot EPI: Parallel Imaging Reconstruction

(A) MUSE on 2-shot EPI


$b = 1000 \text{ s/mm}^2$



slice #28




slice #39




slice #48

$b = 2500 \text{ s/mm}^2$



## 2-Shot EPI: Overlapping LLR Regularization

(B) Our joint recon on 2-shot EPI



## Conclusion

## Conclusion

- ▶ Undersampled iEPI acquisition at ultra-high field

## Conclusion

- ▶ Undersampled iEPI acquisition at ultra-high field
- ▶ Joint  $k$ - $q$ -space reconstruction
- ▶ Spatial-diffusion overlapping LLR regularization

## Conclusion

- ▶ Undersampled iEPI acquisition at ultra-high field
- ▶ Joint  $k$ - $q$ -space reconstruction
- ▶ Spatial-diffusion overlapping LLR regularization
- ▶ Preliminary in vivo results show better image quality compared to MUSE

## Conclusion

- ▶ Undersampled iEPI acquisition at ultra-high field
- ▶ Joint  $k$ - $q$ -space reconstruction
- ▶ Spatial-diffusion overlapping LLR regularization
- ▶ Preliminary in vivo results show better image quality compared to MUSE
- ▶ Please refer to:  
[https://github.com/ZhengguoTan/demo\\_jets\\_diffusion\\_mri\\_7t.git](https://github.com/ZhengguoTan/demo_jets_diffusion_mri_7t.git)  
for future updates!

## Conclusion

- ▶ Undersampled iEPI acquisition at ultra-high field
- ▶ Joint  $k$ - $q$ -space reconstruction
- ▶ Spatial-diffusion overlapping LLR regularization
- ▶ Preliminary in vivo results show better image quality compared to MUSE
- ▶ Please refer to:  
[https://github.com/ZhengguoTan/demo\\_jets\\_diffusion\\_mri\\_7t.git](https://github.com/ZhengguoTan/demo_jets_diffusion_mri_7t.git)  
for future updates!
- ▶ MR-Pub III at Reproducible Research Study Group at 13:45 on Thursday

## Conclusion

- ▶ Undersampled iEPI acquisition at ultra-high field
- ▶ Joint  $k$ - $q$ -space reconstruction
- ▶ Spatial-diffusion overlapping LLR regularization
- ▶ Preliminary in vivo results show better image quality compared to MUSE
- ▶ Please refer to:  
[https://github.com/ZhengguoTan/demo\\_jets\\_diffusion\\_mri\\_7t.git](https://github.com/ZhengguoTan/demo_jets_diffusion_mri_7t.git)  
for future updates!
- ▶ MR-Pub III at Reproducible Research Study Group at 13:45 on Thursday
- ▶ Thank You for Your Attention!